

## Summary

Vendor neutral data processing and reporting software for oligonucleotide impurity analysis

Assessment of novel progressive deconvolution data processing approach

Different relative quantitation approaches compared including UV, MS, and UV-MS

Accelerating oligonucleotide LC-MS data analysis and reporting

## Introduction

Reliable detection and relative quantitation of synthetic oligonucleotide impurities by LC-MS is essential for process control. The data processing approach used directly impacts impurity detection coverage and relative quantitation consistency. Here we compare two approaches, **Trace Peak Centric (TPC)** and **Progressive Deconvolution (PD)**, using Protein Metrics Byos Oligo software.

Figure 1: Trace Peak Centric (TPC) Data Processing Approach





Figure 2: Progressive Deconvolution (PD) Data Processing Approach




## Methods

Triplicate Agilent DNA 40mer LC-HRMS datafiles were analyzed using **Byos Oligo (v6.0)** with TPC and PD data processing workflows (Figures 1 and 2). Monoisotopic masses were matched to the FLP and defined impurity types (clips, shortmers, longmers, substitutions, and chemical modifications) within 10 ppm. Relative quantitation was estimated using UV-only, UV-MS, and MS-only approaches. Automated report generation was performed using Byos template-driven reporting.

## Results – UV relative quantitation (TPC)

Figure 3: UV Trace Peak Integration and % Relative Peak Area




10x UV trace peaks were found per sample and integrated using Byos AutoCompute trace peak finder.

Byos automatically calculates relative UV peak area, estimating 90% FLP purity (RSD < 10%).

## Results – UV-MS relative quantitation (TPC)

Figure 4: UV-MS calculation example



To account for coeluting impurities, Byos automatically adjusts relative UV peak area based on MS intensity ratios calculated for each trace peak (UV-MS). Using this approach, we estimate 70% FLP purity (RSD < 10%).

Figure 5: UV-MS FLP and top 5 impurities

Oligo name	Sample				
	Replicate	1	2	3	[Avg]   [RSD]
Agilent 40-mer Reference		70.62	67.46	71.63	69.90   3.11
Agilent 40-mer n-1 Ade		2.85	3.17	2.62	2.88   9.58
Agilent 40-mer n-2 Ade		2.47	2.67	2.46	2.53   4.61
Agilent 40-mer n-3 Ade		2.20	2.42	2.00	2.21   9.65
Agilent 40-mer 11-40(A,dG,dT,dG,dA,...)		1.58	1.83	1.54	1.65   9.67
Agilent 40-mer 12-40(A,dG,dT,dG,dA,...)		1.55	1.79	1.53	1.62   8.87

## Results – MS relative quantitation (TPC and PD)


Relative quantitation using deconvoluted mass Intensity (TPC) and XIC AUC (PD). Intensity yielded superior reproducibility (%RSD < 20%). The lower reproducibility observed with PD XIC AUC is likely due to limited MS scans per chromatographic peak; faster acquisition rates or broader LC peaks would improve reliability.

Figure 6: MS only relative quantitation TPC vs. PD (top 10 masses detected by abundance)

Oligonucleotide	Mono mass	Trace peak centric					Progressive deconvolution				
		Rep 1	Rep 2	Rep 3	Avg	%RSD	Rep 1	Rep 2	Rep 3	Avg	%RSD
Agilent 40-mer FLP	12286.11	39.21%	37.20%	39.79%	38.73%	3.51%	38.38%	36.42%	37.35%	37.38%	2.62%
Clip (34-40)	2079.40	3.31%	3.33%	3.45%	3.36%	2.13%	3.79%	3.89%	5.55%	4.41%	22.40%
Clip (34-40)-3'PO	2199.36	2.54%	2.62%	2.41%	2.52%	4.29%	3.43%	3.44%	3.28%	3.39%	2.69%
Clip (29-40)	3643.65	2.19%	2.39%	2.15%	2.24%	5.76%	2.45%	2.52%	2.34%	2.44%	3.87%
Clip (11-40)	9259.99	2.17%	2.25%	2.17%	2.20%	2.09%	0.71%	0.78%	1.23%	0.90%	31.56%
Clip (17-40)	7382.28	2.00%	2.15%	1.99%	2.04%	4.27%	2.39%	2.55%	2.36%	2.43%	4.02%
Clip (25-40)	4878.85	1.93%	1.95%	1.99%	1.96%	1.57%	1.68%	1.71%	1.06%	1.48%	24.68%
Clip (21-40)	6138.07	1.98%	1.99%	1.85%	1.94%	3.86%	1.45%	1.47%	1.30%	1.41%	6.57%
Clip (25-40)-3'PO	4958.82	1.76%	1.80%	1.80%	1.79%	1.07%	1.99%	1.97%	1.94%	1.97%	1.30%

## Results – Comparing TPC and PD

Figure 7: Impurities detected TPC vs. PD



TPC and PD: Detect FLP and >100 impurity masses (below 0.01% Rel MS abundance), with rel. abundance estimated for each identified mass.

**TPC Benefits:** Enables UV, MS, and UV-MS based quantitation. UV-MS provides a hybrid approach with advantages over UV-only and MS-only methods. Simpler, faster processing makes TPC well-suited for high-throughput analysis.

**PD Benefits:** Signal processing limits noise, enabling detection of very low abundance impurities. A single mass assignment per impurity simplifies interpretation and avoids duplicate assignments seen across multiple trace peaks.

## Conclusion

TPC and PD approaches both successfully identified the FLP and >100 impurity masses, including n-1, n+1, clips (6-39mer), and other modifications. UV-only, MS-only, and UV-MS quantitation provide complementary approaches for estimating relative abundance. The analysis and reporting are automated using Byos Oligo. For this dataset PD demonstrates higher sensitivity (+31 impurity mass assignments), while TPC relative quantitation showed superior reproducibility compared to PD, which would benefit from more MS scans per chromatographic peak. TPC and PD offer unique and complementary advantages.