

Summary

Integrated Framework for System Suitability

Combines automated dashboards and advanced visualization tools to ensure LC-MS instrument readiness and data confidence.

Automation Enhances Efficiency and Compliance

Real-time monitoring via Protein Metrics' Byosphere platform reduces manual intervention, accelerates troubleshooting, and supports regulatory data integrity.

Proactive Validation Improves Data Quality

Early detection of anomalies mitigates risks in biopharmaceutical characterization, ensuring reproducibility and robust decision-making in development pipelines.


Introduction

Reliable LC-MS data generation for biopharmaceutical characterization requires rigorous system suitability assessment before analysis. This work introduces an integrated framework using automated dashboards and visualization tools to ensure instrument readiness and data confidence. Suitability protocols monitor critical indicators such as retention time stability, mass accuracy, peptide digest quality, and chemical quality attributes (CQAs) like oxidation and deamidation, enabling early detection of deviations. Leveraging Protein Metrics' Byosphere platform, automated processing and dashboarding provide real-time evaluation across instruments and workflows. Browser-based access supports global review, while Deep Query dashboards streamline anomaly investigation and root cause analysis, reducing downtime and manual intervention.

This approach improves reproducibility, mitigates risks from inaccurate quantitation or incomplete peptide coverage, and supports regulatory compliance through automated reporting. Case studies show enhanced efficiency and data quality via early anomaly detection and corrective action. Routine implementation strengthens confidence in results, enabling robust biopharmaceutical characterization and informed decision-making. Combining system suitability with automation and visualization ensures high-quality, reproducible LC-MS data in modern analytical environments.

<ul style="list-style-type: none"> Molecular weight confirmation Clipping Covalent modifications 	<ul style="list-style-type: none"> Molecular weight confirmation Glycoform characterization Clipping PTM characterization 	<ul style="list-style-type: none"> Molecular weight confirmation Glycation Clipping Proline amidation Covalent modifications 	<ul style="list-style-type: none"> Aggregation Noncovalent interactions Excess chain expression 	<ul style="list-style-type: none"> Trypsin and Chymotrypsin digests Sequence variant analysis PTM characterization Sequence confirmation Contamination
Intact Mass	Reduced Mass	Deglycosylated Intact Mass	Native SEC-MS	Peptide Mapping




Automated Cloud-Based Data Processing & Analysis



Globally Accessible & Notifications

The screenshot shows the 'Byosphere: Virtual Client' interface. It features a dashboard with various system suitability metrics and a 'Background Alerts' panel. The alerts panel includes a 'Mass Error Alert' with a description of the error (e.g., 'Mass error > 3 PPM') and an 'Alert Filter' section. The interface is designed for global accessibility and provides email notifications for system suitability dashboards.

Evaluate System Suitability Controls Earlier-Automation allows us to shift the paradigm



Beyond Quality Control

